

Product datasheet for **HP300330**

MIR365-2 Human qPCR Primer Pair (MI0000769)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100126356
Forward Sequence:	GGGACTTTCAGGGGCA
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0009199
Synonyms:	hsa-miR-365; MI0000769; MIMAT0009199; MIR365-2 hsa-mir-365-2
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20 °C.



[View online »](#)