

Product datasheet for **HP300692**

MIR1269b Human qPCR Primer Pair (MI0016888)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Locus ID: | 100616494 |
| Forward Sequence: | CTGGACTGAGCCATGCTACT |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0019059 |
| Synonyms: | hsa-miR-1269b; MI0016888; MIMAT0019059; MIR1269b hsa-mir-1269b |
| Components: | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | The primer mix is stable for one year from date of shipping. Store at -20°C. |



[View online »](#)