

## Product datasheet for **HP300787**

### MIR3179-1 Human qPCR Primer Pair (MI0014213)

#### Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100422960
Forward Sequence:	AGAAGGGTGAAATTTAAACG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0015056
Synonyms:	hsa-miR-3179; MI0014213; MIMAT0015056; MIR3179-1 hsa-mir-3179-1
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)