

## Product datasheet for **HP300946**

### MIR4421 Human qPCR Primer Pair (MI0016758)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Locus ID:         | 100616189  |
| Forward Sequence: | ACCTGTCTGTGGAAGGAGC  |
| Reverse Sequence: | GAACATGTCTGCGTATCTC  |
| ACCN:             | MIMAT0018934   |
| Synonyms:         | hsa-miR-4421; MI0016758; MIMAT0018934; MIR4421 hsa-mir-4421  |
| Components:       | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | The primer mix is stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)