

Product datasheet for **HP300964**

MIR4446 Human qPCR Primer Pair (MI0016789)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Locus ID: | 100616476 |
| Forward Sequence: | ATTTCCTGCCATTCCCTTG |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0019233 |
| Synonyms: | hsa-miR-4446-5p; MI0016789; MIMAT0019233; MIR4446 hsa-mir-4446 |
| Components: | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | The primer mix is stable for one year from date of shipping. Store at -20°C. |



[View online »](#)