

Product datasheet for **HP301047**

MIR4653 Human qPCR Primer Pair (MI0017281)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100616117
Forward Sequence:	TCTCTGAGCAAGGCTTAAC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0019718
Synonyms:	hsa-miR-4653-5p; MI0017281; MIMAT0019718; MIR4653 hsa-mir-4653
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)