

## Product datasheet for **HP301225**

### MIR550a-3 Human qPCR Primer Pair (MI0003762)

#### Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100616354
Forward Sequence:	AGTGCCTGAGGGAGTAAG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0020925
Synonyms:	hsa-miR-550a-3-5p; MI0003762; MIMAT0020925; MIR550a-3 hsa-mir-550a-3
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)