

## Product datasheet for **HP301234**

### **MIR5588 Human qPCR Primer Pair (MI0019147)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Forward Sequence:</b>	ACTGGCATTAGTGGGACTT
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0022295
<b>Synonyms:</b>	hsa-miR-5588-5p; MI0019147; MIMAT0022295; MIR5588 hsa-mir-5588
<b>Components:</b>	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)