

## Product datasheet for **HP301271**

### **MIR6080 Human qPCR Primer Pair (MI0020357)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Forward Sequence:</b>	TCTAGTGCGGCGTTCC
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0023705
<b>Synonyms:</b>	hsa-miR-6080; MI0020357; MIMAT0023705; MIR6080 hsa-mir-6080
<b>Components:</b>	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)