

## Product datasheet for HPP6009D

### Primer Panel

#### Product data:

Product Type: qPCR Primer Panels  
Description: T-Cell/B-Cell Activation (Human) qPCR primer panel - 96 X 200 rxns in matrix tubes  
panel: T-Cell/B-Cell Activation  
panel\_link: /assets/documents/qpcr/t-cellb-cell.xls

#### Product images:

	1	2	3	4	5	6	7	8	9	10	11	12
A	1	2	3	4	5	6	7	8	9	10	11	12
B	13	14	15	16	17	18	19	20	21	22	23	24
C	25	26	27	28	29	30	31	32	33	34	35	36
D	37	38	39	40	41	42	43	44	45	46	47	48
E	49	50	51	52	53	54	55	56	57	58	59	60
F	61	62	63	64	65	66	67	68	69	70	71	72
G	73	74	75	76	77	78	79	80	81	82	83	84
H	85	86	87	88	89	90	91	92	93	94	95	96



[View online »](#)