

Product datasheet for HPP6017D

Primer Panel

Product data:

Product Type: qPCR Primer Panels
Description: Signal Transduction (Human) qPCR primer panel - 96 X 200 rxns in matrix tubes
panel: Signal Transduction
panel_link: /assets/documents/qpcr/signaltransduction.xls

Product images:

	1	2	3	4	5	6	7	8	9	10	11	12
A	1	2	3	4	5	6	7	8	9	10	11	12
B	13	14	15	16	17	18	19	20	21	22	23	24
C	25	26	27	28	29	30	31	32	33	34	35	36
D	37	38	39	40	41	42	43	44	45	46	47	48
E	49	50	51	52	53	54	55	56	57	58	59	60
F	61	62	63	64	65	66	67	68	69	70	71	72
G	73	74	75	76	77	78	79	80	81	82	83	84
H	85	86	87	88	89	90	91	92	93	94	95	96



[View online »](#)