

Product datasheet for **HPRT301**

TissueScan, Prostate Cancer cDNA Array I

Product data:

Product Type:	TissueScan
Array Category:	Cancer Specific
Clinical Info:	/assets/documents/tissuescan/hprt301.xls
Diagnosis Category:	Prostate Cancer
Tissue:	Prostate
Array Name:	Prostate Cancer cDNA Array I


Product images:




TissueScan product image.



[View online »](#)



A Human Prostate Cancer real time PCR Gene Expression panel used in real-time PCR. The beta-actin contents in the 48 human cDNAs were detected using SYBR Green I protocol and performed in PE7700.
Fluorescence reading - cycle number



A Human Prostate Cancer real time PCR Gene Expression panel was used in SYBR Green I based real-time PCR. The expression profile of the human Hepsin transcript (ACCN#: NM_182983) was detected using SYBR Green I protocol and performed in PE7700.