

## Product datasheet for **HPRT502**

### TissueScan, Prostate Cancer cDNA Array II

#### Product data:

Product Type:	TissueScan
Array Category:	Cancer Specific
Clinical Info:	/assets/documents/tissuescan/hprt502.xls
Diagnosis Category:	Prostate Cancer
Tissue:	Prostate
Array Name:	Prostate Cancer cDNA Array II


#### Product images:



TissueScan product image.




[View online »](#)



A Human Prostate Cancer real time PCR Gene Expression panel used in real-time PCR. The beta-actin contents in the 48 human cDNAs were detected using SYBR Green I protocol and performed in PE7700.

Fluorescence reading - cycle number.



A Human Prostate Cancer real time PCR Gene Expression panel was used in SYBR Green I based real-time PCR. The expression profile of the human Hepsin (NM\_182983) transcript was detected using SYBR Green I protocol and performed in PE7700.

Normalized relative expression among 48 Prostate samples