

Product datasheet for **HSRT101**

TissueScan, Sarcoma cDNA Array I

Product data:

Product Type:	TissueScan
Array Category:	Cancer Specific
Clinical Info:	/assets/documents/tissuescan/hsrt101.xls
Diagnosis:	Sarcoma Cancer
Tissue:	Soft tissue, Muscle, Bone, Cartilage
Array Name:	Sarcoma cDNA Array I


Product images:



TissueScan product image.




[View online »](#)



A human sarcoma real time PCR Gene Expression panel used in real time PCR. The beta-actin contents in the 48 human cDNAs were detected using SYBR Green I protocol and performed in a PE7700 .

Fluorescence reading - cycle number



A human sarcoma real time PCR Gene Expression panel used in SYBR Green based real time PCR. The expression profile of the human Survivin (NM_001168) was detected using SYBR Green I protocol and performed in a PE7700 .

Normalized relative expression among 48 sarcoma samples