

Product datasheet for **KN200005RB**

DLL3 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	DLL3
Locus ID:	10683
Components:	KN200005G1 , DLL3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN200005G2 , DLL3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN200005RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_016941 , NM_203486
UniProt ID:	Q9NYJ7
Synonyms:	SCDO1
Summary:	This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Mutations in this gene cause autosomal recessive spondylocostal dysostosis 1. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

