

## Product datasheet for **KN200364LP**

### FCGRT Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	FCGRT
Locus ID:	2217
Components:	<b>KN200364G1</b> , FCGRT gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN200364G2</b> , FCGRT gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN200364LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001136019</a> , <a href="#">NM_004107</a>
UniProt ID:	<a href="#">P55899</a>
Synonyms:	alpha-chain; FCRN
Summary:	This gene encodes a receptor that binds the Fc region of monomeric immunoglobulin G. The encoded protein transfers immunoglobulin G antibodies from mother to fetus across the placenta. This protein also binds immunoglobulin G to protect the antibody from degradation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009]



[View online »](#)

## Product images:

