

Product datasheet for **KN200714LP**

ICAM1 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	ICAM1
Locus ID:	3383
Components:	KN200714G1 , ICAM1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN200714G2 , ICAM1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN200714LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_000201
UniProt ID:	P05362
Synonyms:	BB2; CD54; P3.58
Summary:	This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

