

Product datasheet for **KN201052LP**

CEP55 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	CEP55
Locus ID:	55165
Components:	KN201052G1 , CEP55 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN201052G2 , CEP55 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN201052LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001127182 , NM_018131
UniProt ID:	Q53EZ4
Synonyms:	C10orf3; CT111; URCC6
Summary:	Plays a role in mitotic exit and cytokinesis (PubMed:16198290, PubMed:17853893). Recruits PDCD6IP and TSG101 to midbody during cytokinesis. Required for successful completion of cytokinesis (PubMed:17853893). Not required for microtubule nucleation (PubMed:16198290). Plays a role in the development of the brain and kidney (PubMed:28264986).[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

