

## Product datasheet for **KN201175RB**

### GIPC (GIPC1) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	GIPC
Locus ID:	10755
Components:	<b>KN201175G1</b> , GIPC gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN201175G2</b> , GIPC gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN201175RBD</b> , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_005716</a> , <a href="#">NM_202467</a> , <a href="#">NM_202468</a> , <a href="#">NM_202469</a> , <a href="#">NM_202470</a> , <a href="#">NM_202494</a>
UniProt ID:	<a href="#">O14908</a>
Synonyms:	C19orf3; GIPC; GLUT1CBP; Hs.6454; IIP-1; NIP; RGS19IP1; SEMCAP; SYNECTIIN; SYNECTIN; TIP-2
Summary:	GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking (Lee et al., 2008 [PubMed 18775991]).[supplied by OMIM, Apr 2009]



[View online »](#)

## Product images:

