

## Product datasheet for **KN201403RB**

### ABCF3 Human Gene Knockout Kit (CRISPR)

#### Product data:

**Product Type:** Knockout Kits (CRISPR)

**Format:** 2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control

**Donor DNA:** RFP-BSD

**Symbol:** ABCF3

**Locus ID:** 55324

**Components:** **KN201403G1**, ABCF3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)  
**KN201403G2**, ABCF3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)  
**KN201403RBD**, donor DNA containing left and right homologous arms and RFP-BSD functional cassette.  
**GE100003**, scramble sequence in pCas-Guide vector

**Disclaimer:** These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

**RefSeq:** [NM\\_018358](#), [NM\\_001351298](#), [NM\\_001351299](#), [NM\\_001351300](#)

**UniProt ID:** [Q9NUQ8](#)

**Synonyms:** EST201864

**Summary:** This gene encodes a member of the ATP-binding cassette (ABC) transporter superfamily. ATP-binding cassette proteins transport various molecules across extra- and intracellular membranes. The protein encoded by this gene displays antiviral effect against flaviviruses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2017]



[View online »](#)

## Product images:

