

## Product datasheet for **KN201616LP**

### **S100A4 Human Gene Knockout Kit (CRISPR)**

#### **Product data:**

<b>Product Type:</b>	Knockout Kits (CRISPR)
<b>Format:</b>	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
<b>Donor DNA:</b>	Luciferase-Puro
<b>Symbol:</b>	S100A4
<b>Locus ID:</b>	6275
<b>Components:</b>	<b>KN201616G1</b> , S100A4 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN201616G2</b> , S100A4 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN201616LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
<b>Disclaimer:</b>	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
<b>RefSeq:</b>	<a href="#">NM_002961</a> , <a href="#">NM_019554</a>
<b>UniProt ID:</b>	<a href="#">P26447</a>
<b>Synonyms:</b>	18A2; 42A; CAPL; FSP1; MTS1; P9KA; PEL98
<b>Summary:</b>	The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in motility, invasion, and tubulin polymerization. Chromosomal rearrangements and altered expression of this gene have been implicated in tumor metastasis. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]



[View online »](#)

## Product images:

