

Product datasheet for **KN201807RB**

Glucose 6 Phosphate Dehydrogenase (G6PD) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Glucose 6 Phosphate Dehydrogenase
Locus ID:	2539
Components:	KN201807G1 , Glucose 6 Phosphate Dehydrogenase gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN201807G2 , Glucose 6 Phosphate Dehydrogenase gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN201807RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_000402](#), [NM_001042351](#), [NM_001360016](#)

UniProt ID: [P11413](#)

Synonyms: G6PD1

Summary: This gene encodes glucose-6-phosphate dehydrogenase. This protein is a cytosolic enzyme encoded by a housekeeping X-linked gene whose main function is to produce NADPH, a key electron donor in the defense against oxidizing agents and in reductive biosynthetic reactions. G6PD is remarkable for its genetic diversity. Many variants of G6PD, mostly produced from missense mutations, have been described with wide ranging levels of enzyme activity and associated clinical symptoms. G6PD deficiency may cause neonatal jaundice, acute hemolysis, or severe chronic non-spherocytic hemolytic anemia. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

