

Product datasheet for **KN201828LP**

SMAD4 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	SMAD4
Locus ID:	4089
Components:	KN201828G1 , SMAD4 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN201828G2 , SMAD4 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN201828LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette.

Homologous arm and Luciferase-Puro sequences:

pUC vector backbone in gray; **Left arm sequence in blue**; **Luciferase-Puro in green**; **Right arm in violet**

```
AAGGCGAGTT ACATGATCCC CCATGTTGTG CAAAAAGCG GTTAGCTCCT TCGGTCCTCC GATCGTTGTC
AGAAGTAAGT TGGCCGAGT GTTATCACTC ATGGTTATGG CAGCACTGCA TAATTCTCTT ACTGTCAATGC
CATCCGTAAG ATGCTTTTCT GTGACTGGTG AGTACTCAAC CAAGTCATTC TGAGAATAGT GTATGCGGCG
ACCGAGTTGC TCTTGCCCGG CGTCAATACG GGATAATACC GCGCCACATA GCAGAACTTT AAAAGTGCTC
ATCATTGGAA AACGTTCTTC GGGGCGAAAA CTCTCAAGGA TCTTACCCTG GTTGAGATCC AGTTCGATGT
AACCCACTCG TGCACCCAAC TGATCTTCAG CATCTTTTAC TTTCACCAGC GTTTCTGGGT GAGCAAAAAAC
AGGAAGGCAA AATGCCGCAA AAAAGGGAAT AAGGGCGACA CGGAAATGTT GAATACTCAT ACTCTTCCTT
TTTCAATATT ATTGAAGCAT TTATCAGGT TATTGTCTCA TGAGCGGATA CATATTTGAA TGTATTTAGA
AAAATAACA AATAGGGGTT CCGCGCAT TCCCCGAAA AGTGCCACCT GACGTCTAAG AAACCATTAT
TATCATGACA TTAACCTATA AAAATAGGCG TATCACGAGG CCTTTTCGGG TCGCGCGTTT CGGTGATGAC
GGTGAAAACC TCTGACACAT GCAGCTCCCG TTGACGGTCA CAGCTTGCTT GTAAGCGGAT GCCGGGAGCA
GACAAGCCCG TCAGGGCGCG TCAGCGGGTG TTGGCGGGTG TCGGGGCTGG CTAACTATG CGGCATCAGA
GCAGATTGTA CTGAGAGTGC ACCATAAAT TGTAACGTT AATATTTTGT TAAAATTCGC GTTAAATTTT
TGTTAAATCA GTCATTTTT TAACCAATAG GCCGAAATCG GCAAAATCCC TTATAATCA AAAGAATAGC
CCGAGATAGG GTTGAGTGTT GTTCCAGTTT GGAACAAGAG TCCACTATTA AAGAACGTGG ACTCCAACGT
CAAAGGGCGA AAAACCGTCT ATCAGGGCGA TGGCCACTA CGTGAACCAT CACCAAATC AAGTTTTTTG
GGTTCGAGGT GCCGTAAGC ACTAAATCGG AACCTAAAG GGAGCCCCG ATTTAGAGT TGACGGGGAA
AGCCGGCGAA CGTGCGGAGA AAGGAAGGGA AGAAAGCGAA AGGAGCGGGC GCTAGGGCGC TGGCAAGTGT
AGCGGTACG CTGCGCGTAA CCACCACACC CGCCGCGCTT AATGCGCCG TACAGGGCGC GACTATGTT
TGCTTTGACG TATGCGGTGT GAAATACCGC ACAGATGCGT AAGGAGAAAA TACCGCATCA GGCGCCATTC
GCCATTACAG CTGCGCAACT GTTGGGAAGG GCGATCGGTG CGGGCCTCTT CGCTATTACG CCAGCTGGCG
AAAGGGGAT GTGCTGCAAG GCGATTAAGT TGGGTAACGC CAGGGTTTTT CCAGTACGCA CGTTGTAAAA
CGACGGCCAG TGAATTGGAG GCTACAGTCA GTGGAGAGGA CTTTCACAGG CTGTCGCCGT GCTCATTTGA
TAACTGCCG TTATTCATGC GACACTTTTT TAATTGGTCT CTTTAAACA AGGTGAGTGT TAAGTCTAG
GTAGAATTC ATACCAGCG ATGGTAAAT TTTCATAAT ATTAGTAATA GACTTTTAGG TTTAACTAGT
```




```

CCTGCTCTTT ACTGAAGGCT CTTTACTATT GCTTTATGAT AATGTTTCAT AGTTGGATAT CATAATTTAA
ACAAGCAAAA CCAAATTAAG GGCCAGCTCA TTCCTCCCAC TCATGATCTA TAGATCTATA GATCTCTCGT
GGGATCATTG TTTTCTCTT GATTCCCCTT TTGTGGTTCT AAGTACTGTG GTTTCCAAAT GTGTCAGTTT
CATAGCCTGA AGAACGAGAT CAGCAGCCTC TGTTCCACAT ACAC TTCATT CTCAGTATTG TTTTGCCAAG
TTCTAATTCC ATCAGAAGCT GGTTCGAGATC CGGAACCCTT AATATAACTT CGTATAATGT ATGCTATACG
AAGTTATTAG GTCCTCGAA GAGGTTCACT AGGCGCGCCA TTTGCAAAA GAGCAATTGA AAGTTTGGTA
AAGAAGCTGA AGGAGAAAA AGATGAATTG GATTCTTTAA TAACAGCTAT AACTACAAAT GGAGCTCATC
CTAGTAAATG TGTTACCATA CAGAGAACAT TGGATGGGAG GCTTCAGGTT AGTCTTATAA GAGTTTTTCT
ATACCCTCTA TGGTGGCAGA TTTAAAAACT TGCTACGTTT CCTTTCAGC TACTACAGGG TAATTTAATT
TGTGCTCCAT CTCTTCAGAT ACTGTGCATC CTGTACAAAT ATGCATTATG GGAATTTCTG GAAGAATATG
CAAAAGATAA AATCTTAAAG TTTTTTAATG TTCTACTCAG AAAAATGTTC AATGGAGAAA ATTTGGAAAA
TAAAAATAAA AAGCAGAAAA CCCACACTTC GCAAGAATA CCTTTAGCAC TTTGGTATTG ATTTTTGTAG
CCTGTTGTTT TCTTTTTTTT TTTTTTCCTA GACGGAGTCT TGCTCTGTCA CCCAGGCTGG AGTACAGTGG
CGTGATCTTG GCCCACTGCA ACCTCCACCT CTGGGGTCA AGCAATCCTC CTGCCTCAGC CTCCTAAGTA
GCTGGGACTT CACTCTCGCC GGTGGACTT TAGATCAGAA GGGATCTTGC TGCCGCCGA AAGAGGAAGG
GCTGGAAGAG GAAGGAGCTT GCGTAATCA TGTCATAGC TGTTTCCTGT GTGAAATGT TATCCGCTCA
CAATTCACA CAACATACGA GCCGGAAGCA TAAAGTGTA AGCCTGGGGT GCCTAATGAG TGAGCTAACT
CACATTAATT GCGTTGCGCT CACTGCCCGC TTTCCAGTCG GAAACCTGT CGTGCCAGCT GCATTAATGA
ATCGGCCAAC GCGCGGGGAG AGGCGGTTTG CGTATTGGGC GCTCTTCCGC TTCCTCGCTC ACTGACTCGC
TGCCTCGGT CGTTCGGCTG CCGCGAGCGG TATCAGCTCA CTCAAAGGCG GTAATACGGT TATCCACAGA
ATCAGGGGAT AACGCAGGAA AGAACATGTG AGCAAAAGC CAGCAAAAGG CCAGGAACCG TAAAAAGGCC
GCGTTGCTGG CGTTTTTCCA TAGGCTCCGC CCCCTGACG AGCATCACAA AAATCGACCG TCAAGTCAGA
GGTGCGGAAA CCCGACAGGA CTATAAAGAT ACCAGGCGTT TCCCCGGA AGCTCCCTCG TGCCTCTCC
TGTCCGACC CTGCCGCTTA CCGGATACCT GTCCGCTTT CTCCCTTCGG GAAGCGTGGC GCTTTCTCAT
AGCTCACGCT GTAGGTATCT CAGTTCGGTG TAGGTCGTTT GCTCCAAGCT GGGCTGTGTG CACGAACCCC
CCGTTACGCC CGACCGCTGC GCCTTATCCG GTAACTATCG TCTTGAGTCC AACCCGGTAA GACACGACTT
ATCGCCACTG GCAGCAGCCA CTGGTAACAG GATTAGCAGA GCGAGGTATG TAGGCGGTGC TACAGAGTTC
TTGAAGTGGT GGCCTAACTA CCGCTACACT AGAAGAACAG TATTTGGTAT CTGCGCTCTG CTGAAGCCAG
TTACCTTCGG AAAAAGAGTT GGTAGCTCTT GATCCGGCAA ACAACCACC GCTGGTAGCG GTGGTTTTTT
TGTTTGCAAG CAGCAGATTA CCGCAGAAA AAAAGGATCT CAAGAAGATC CTTTGATCTT TTCTACGGGG
TCTGACGCTC AGTGGAACGA AACTCACGT TAAGGGATTT TGTCATGAG ATTATCAAAA AGGATCTTCA
CCTAGATCCT TTTAAATTA AAATGAAGTT TTAATCAAT CTAAAGTATA TATGAGTAAA CTTGGTCTGA
CAGTTACCAA TGCTTAATCA GTGAGGCACC TATCTCAGCG ATCTGTCTAT TTCGTTTCATC CATAGTTGCC
TGAATCCCG TCGTGTAGAT AACTACGATA CGGGAGGGCT TACCATCTGG CCCAGTGTG GCAATGATAC
CGCGAGAACC ACGCTACCG GCTCCAGATT TATCAGCAAT AAACCAGCCA GCCGGAAGGG CCGAGCGCAG
AAGTGGTCTT GCAACTTTAT CCGCTCCAT CCAGTCTATT AATTGTTGCC GGAAGCTAG AGTAAGTAGT
TCGCCAGTTA ATAGTTTGCG CAACGTTGTT GCCATTGCTA CAGGCATCGT GGTGTACCGC TCGTCGTTT
GTATGGCTTC ATTCAGCTCC GGTCCCAAC GATC

```

GE100003, scramble sequence in pCas-Guide vector

Disclaimer:

These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq:

[NM_005359](#)

UniProt ID:

[Q13485](#)

Synonyms:

DPC4; JIP; MADH4; MYHRS

Summary:

This gene encodes a member of the Smad family of signal transduction proteins. Smad proteins are phosphorylated and activated by transmembrane serine-threonine receptor kinases in response to transforming growth factor (TGF)-beta signaling. The product of this gene forms homomeric complexes and heteromeric complexes with other activated Smad proteins, which then accumulate in the nucleus and regulate the transcription of target genes. This protein binds to DNA and recognizes an 8-bp palindromic sequence (GTCTAGAC) called the Smad-binding element (SBE). The protein acts as a tumor suppressor and inhibits epithelial cell proliferation. It may also have an inhibitory effect on tumors by reducing angiogenesis and increasing blood vessel hyperpermeability. The encoded protein is a crucial component of the bone morphogenetic protein signaling pathway. The Smad proteins are subject to complex regulation by post-translational modifications. Mutations or deletions in this gene have been shown to result in pancreatic cancer, juvenile polyposis syndrome, and hereditary hemorrhagic telangiectasia syndrome. [provided by RefSeq, Aug 2017]

Product images:
