

## Product datasheet for KN202000RB

**CD9 Human Gene Knockout Kit (CRISPR)** 

**Product data:** 

**Product Type:** Knockout Kits (CRISPR)

Format: 2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control

**Donor DNA:** 

Symbol: CD9

928 Locus ID:

**KN202000G1**, CD9 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) Components:

**KN202000G2**, CD9 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)

KN202000RBD, donor DNA containing left and right homologous arms and RFP-BSD

functional cassette.

GE100003, scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is

> designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the

experimental process.

RefSeq: NM 001769, NM 001330312

**UniProt ID:** P21926

Synonyms: BTCC-1; DRAP-27; MIC3; MRP-1; TSPAN-29; TSPAN29

Summary: This gene encodes a member of the transmembrane 4 superfamily, also known as the

> tetraspanin family. Tetraspanins are cell surface glycoproteins with four transmembrane domains that form multimeric complexes with other cell surface proteins. The encoded protein functions in many cellular processes including differentiation, adhesion, and signal transduction, and expression of this gene plays a critical role in the suppression of cancer cell

motility and metastasis. [provided by RefSeq, Jan 2011]



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



# **Product images:**

### Donor Vector Edited Chromosome



RFP, Luc, and mBFP will be under native gene promoter