

Product datasheet for **KN202015BN**

ETS1 associated protein II (TDP2) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	ETS1 associated protein II
Locus ID:	51567
Components:	KN202015G1 , ETS1 associated protein II gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN202015G2 , ETS1 associated protein II gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN202015BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_016614](#)

UniProt ID: [O95551](#)

Synonyms: AD022; dj30M3.3; EAP2; EAPII; hTDP2; TTRAP

Summary: This gene encodes a member of a superfamily of divalent cation-dependent phosphodiesterases. The encoded protein associates with CD40, tumor necrosis factor (TNF) receptor-75 and TNF receptor associated factors (TRAFs), and inhibits nuclear factor-kappa-B activation. This protein has sequence and structural similarities with APE1 endonuclease, which is involved in both DNA repair and the activation of transcription factors. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

