

## Product datasheet for **KN202532LP**

### Mesothelin (MSLN) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	Mesothelin
Locus ID:	10232
Components:	<b>KN202532G1</b> , Mesothelin gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN202532G2</b> , Mesothelin gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN202532LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001177355</a> , <a href="#">NM_005823</a> , <a href="#">NM_013404</a>
UniProt ID:	<a href="#">Q13421</a>
Synonyms:	MPF; SMRP
Summary:	This gene encodes a preproprotein that is proteolytically processed to generate two protein products, megakaryocyte potentiating factor and mesothelin. Megakaryocyte potentiating factor functions as a cytokine that can stimulate colony formation of bone marrow megakaryocytes. Mesothelin is a glycosylphosphatidylinositol-anchored cell-surface protein that may function as a cell adhesion protein. This protein is overexpressed in epithelial mesotheliomas, ovarian cancers and in specific squamous cell carcinomas. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed. [provided by RefSeq, Feb 2016]



[View online »](#)

## Product images:

