

Product datasheet for **KN203394BN**

CAB39L Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	CAB39L
Locus ID:	81617
Components:	KN203394G1 , CAB39L gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN203394G2 , CAB39L gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN203394BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001031724 , NM_001079670 , NM_001287337 , NM_001287338 , NM_001287339 , NM_030925
UniProt ID:	Q9H9S4
Synonyms:	antigen MLAA-34; bA103J18.3; calcium binding protein 39-like; FLJ12577; MO2L; MO2L, FLJ12577, bA103J18.3, RP11-103J18.3; MO25-BETA; OTTHUMP00000018410; RP11-103J18.3; sarcoma antigen NY-SAR-79; U937-associated antigen
Summary:	Component of a complex that binds and activates STK11/LKB1. In the complex, required to stabilize the interaction between CAB39/MO25 (CAB39/MO25alpha or CAB39L/MO25beta) and STK11/LKB1 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

