

Product datasheet for **KN203818RB**

Integrin beta 1 (ITGB1) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Integrin beta 1
Locus ID:	3688
Components:	KN203818G1 , Integrin beta 1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN203818G2 , Integrin beta 1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN203818RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_002211](#), [NM_033666](#), [NM_033667](#), [NM_033668](#), [NM_033669](#), [NM_133376](#)

UniProt ID: [P05556](#)

Synonyms: CD29; FNRB; GPIIA; MDF2; MSK12; VLA-BETA; VLAB

Summary: Integrins are heterodimeric proteins made up of alpha and beta subunits. At least 18 alpha and 8 beta subunits have been described in mammals. Integrin family members are membrane receptors involved in cell adhesion and recognition in a variety of processes including embryogenesis, hemostasis, tissue repair, immune response and metastatic diffusion of tumor cells. This gene encodes a beta subunit. Multiple alternatively spliced transcript variants which encode different protein isoforms have been found for this gene. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

