

Product datasheet for **KN204105**

ThPok (ZBTB7B) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 GFP-puro donor, 1 scramble control
Donor DNA:	GFP-puro
Symbol:	ThPok
Locus ID:	51043
Components:	KN204105G1 , ThPok gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN204105G2 , ThPok gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN204105D , donor DNA containing left and right homologous arms and GFP-puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001252406 , NM_001256455 , NM_015872 , NR_045515 , NR_046206 , NR_049765
UniProt ID:	O15156
Synonyms:	c-KROX; CKROX; hcKROX; THPOK; ZBTB15; ZFP-67; ZFP67; ZNF857B
Summary:	This gene encodes a zinc finger-containing transcription factor that acts as a key regulator of lineage commitment of immature T-cell precursors. It is necessary and sufficient for commitment of CD4 lineage, while its absence causes CD8 commitment. It also functions as a transcriptional repressor of type I collagen genes. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jan 2012]



[View online »](#)

Product images:

