

## Product datasheet for **KN204454BN**

### TAK1 (MAP3K7) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	TAK1
Locus ID:	6885
Components:	<b>KN204454G1</b> , TAK1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN204454G2</b> , TAK1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN204454BND</b> , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_003188</a> , <a href="#">NM_145331</a> , <a href="#">NM_145332</a> , <a href="#">NM_145333</a>
UniProt ID:	<a href="#">O43318</a>
Synonyms:	MEKK7; TAK1; TGF1a
Summary:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. This kinase mediates the signaling transduction induced by TGF beta and morphogenetic protein (BMP), and controls a variety of cell functions including transcription regulation and apoptosis. In response to IL-1, this protein forms a kinase complex including TRAF6, MAP3K7P1/TAB1 and MAP3K7P2/TAB2; this complex is required for the activation of nuclear factor kappa B. This kinase can also activate MAPK8/JNK, MAP2K4/MKK4, and thus plays a role in the cell response to environmental stresses. Four alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]



[View online »](#)

## Product images:

