

Product datasheet for **KN204497RB**

HSD3B1 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	HSD3B1
Locus ID:	3283
Components:	KN204497G1 , HSD3B1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN204497G2 , HSD3B1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN204497RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_000862](#), [NM_001328615](#)

UniProt ID: [P14060](#)

Synonyms: 3BETAHSD; HSD3B; HSDB3; HSDB3A; I; SDR11E1

Summary: The protein encoded by this gene is an enzyme that catalyzes the oxidative conversion of delta-5-3-beta-hydroxysteroid precursors into delta-4-ketosteroids, which leads to the production of all classes of steroid hormones. The encoded protein also catalyzes the interconversion of 3-beta-hydroxy- and 3-keto-5-alpha-androstane steroids. [provided by RefSeq, Jun 2016]



[View online »](#)

Product images:

