

Product datasheet for **KN205162BN**

Angiotensin II Type 1 Receptor (AGTR1) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Angiotensin II Type 1 Receptor
Locus ID:	185
Components:	KN205162G1 , Angiotensin II Type 1 Receptor gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN205162G2 , Angiotensin II Type 1 Receptor gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN205162BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_000685](#), [NM_004835](#), [NM_009585](#), [NM_031850](#), [NM_032049](#)

UniProt ID: [P30556](#)

Synonyms: AG2S; AGTR1B; AT1; AT1AR; AT1B; AT1BR; AT1R; AT2R1; HAT1R

Summary: Angiotensin II is a potent vasopressor hormone and a primary regulator of aldosterone secretion. It is an important effector controlling blood pressure and volume in the cardiovascular system. It acts through at least two types of receptors. This gene encodes the type 1 receptor which is thought to mediate the major cardiovascular effects of angiotensin II. This gene may play a role in the generation of reperfusion arrhythmias following restoration of blood flow to ischemic or infarcted myocardium. It was previously thought that a related gene, denoted as AGTR1B, existed; however, it is now believed that there is only one type 1 receptor gene in humans. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Aug 2020]



[View online »](#)

Product images:

