

Product datasheet for **KN205282RB**

EDIL3 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	EDIL3
Locus ID:	10085
Components:	KN205282G1 , EDIL3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN205282G2 , EDIL3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN205282RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001278642 , NM_005711
UniProt ID:	O43854
Synonyms:	DEL1
Summary:	The protein encoded by this gene is an integrin ligand. It plays an important role in mediating angiogenesis and may be important in vessel wall remodeling and development. It also influences endothelial cell behavior. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

