

Product datasheet for **KN205306BN**

MTM1 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	MTM1
Locus ID:	4534
Components:	KN205306G1 , MTM1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN205306G2 , MTM1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN205306BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_000252
UniProt ID:	Q13496
Synonyms:	CNM; MTMX; XLMTM
Summary:	This gene encodes a dual-specificity phosphatase that acts on both phosphotyrosine and phosphoserine. It is required for muscle cell differentiation and mutations in this gene have been identified as being responsible for X-linked myotubular myopathy. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

