

## Product datasheet for **KN205317RB**

### PI 3 Kinase Class 3 (PIK3C3) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	PI 3 Kinase Class 3
Locus ID:	5289
Components:	<b>KN205317G1</b> , PI 3 Kinase Class 3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN205317G2</b> , PI 3 Kinase Class 3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN205317RBD</b> , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001308020</a> , <a href="#">NM_002647</a>
UniProt ID:	<a href="#">Q8NEB9</a>
Synonyms:	hVps34; VPS34; Vps34
Summary:	Catalytic subunit of the PI3K complex that mediates formation of phosphatidylinositol 3-phosphate; different complex forms are believed to play a role in multiple membrane trafficking pathways: PI3KC3-C1 is involved in initiation of autophagosomes and PI3KC3-C2 in maturation of autophagosomes and endocytosis. Involved in regulation of degradative endocytic trafficking and required for the abscission step in cytokinesis, probably in the context of PI3KC3-C2 (PubMed:20643123, PubMed:20208530). Involved in the transport of lysosomal enzyme precursors to lysosomes. Required for transport from early to late endosomes (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

## Product images:

