

Product datasheet for **KN206071RB**

RO60 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	RO60
Locus ID:	6738
Components:	KN206071G1 , RO60 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN206071G2 , RO60 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN206071RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001042369 , NM_001042370 , NM_001173524 , NM_001173525 , NM_004600 , NR_033393 , NM_001331020
UniProt ID:	P10155
Synonyms:	RO60; RORN; SSA2
Summary:	RNA-binding protein that binds to misfolded non-coding RNAs, pre-5S rRNA, and several small cytoplasmic RNA molecules known as Y RNAs. May stabilize some of these RNAs and protect them from degradation (PubMed:18056422). Binds to endogenous Alu retroelements which are induced by type I interferon and stimulate porinflammatoiry cytokine secretion. Regulates the expression of Alu retroelements as well as inflammatory genes (PubMed:26382853).[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

