

Product datasheet for **KN206187BN**

TXNDC (TMX1) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	TXNDC
Locus ID:	81542
Components:	KN206187G1 , TXNDC gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN206187G2 , TXNDC gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN206187BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_030755
UniProt ID:	Q9H3N1
Synonyms:	PDIA11; TMX; TXNDC; TXNDC1
Summary:	This gene encodes a member of the disulfide isomerase (PDI) family of endoplasmic reticulum (ER) proteins that catalyze protein folding and thiol-disulfide interchange reactions. The encoded protein has an N-terminal ER-signal sequence, a catalytically active thioredoxin domain, and one transmembrane domain. Unlike most members of this gene family, it lacks a C-terminal ER-retention sequence. The mature membrane-bound protein can both oxidize and reduce disulfide bonds and acts selectively on membrane-associated polypeptides. [provided by RefSeq, Jan 2017]



[View online »](#)

Product images:

