

## Product datasheet for **KN206431BN**

### AXL Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	AXL
Locus ID:	558
Components:	<b>KN206431G1</b> , AXL gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN206431G2</b> , AXL gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN206431BND</b> , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001278599</a> , <a href="#">NM_001699</a> , <a href="#">NM_021913</a>
UniProt ID:	<a href="#">P30530</a>
Synonyms:	ARK; JTK11; Tyro7; UFO
Summary:	The protein encoded by this gene is a member of the Tyro3-Axl-Mer (TAM) receptor tyrosine kinase subfamily. The encoded protein possesses an extracellular domain which is composed of two immunoglobulin-like motifs at the N-terminal, followed by two fibronectin type-III motifs. It transduces signals from the extracellular matrix into the cytoplasm by binding to the vitamin K-dependent protein growth arrest-specific 6 (Gas6). This gene may be involved in several cellular functions including growth, migration, aggregation and anti-inflammation in multiple cell types. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, Jul 2013]



[View online »](#)

## Product images:

