

## Product datasheet for **KN206453RB**

### CD4 Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	CD4
Locus ID:	920
Components:	<b>KN206453G1</b> , CD4 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN206453G2</b> , CD4 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN206453RBD</b> , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_000616</a> , <a href="#">NM_001195014</a> , <a href="#">NM_001195015</a> , <a href="#">NM_001195016</a> , <a href="#">NM_001195017</a> , <a href="#">NR_036545</a>
UniProt ID:	<a href="#">P01730</a>
Synonyms:	CD4mut
Summary:	This gene encodes the CD4 membrane glycoprotein of T lymphocytes. The CD4 antigen acts as a coreceptor with the T-cell receptor on the T lymphocyte to recognize antigens displayed by an antigen presenting cell in the context of class II MHC molecules. The CD4 antigen is also a primary receptor for entry of the human immunodeficiency virus through interactions with the HIV Env gp120 subunit. This gene is expressed not only in T lymphocytes, but also in B cells, macrophages, granulocytes, as well as in various regions of the brain. The protein functions to initiate or augment the early phase of T-cell activation, and may function as an important mediator of indirect neuronal damage in infectious and immune-mediated diseases of the central nervous system. Multiple alternatively spliced transcript variants encoding different isoforms have been identified in this gene. [provided by RefSeq, May 2020]



[View online »](#)

Product images:

