

Product datasheet for **KN206468BN**

Occludin (OCLN) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Occludin
Locus ID:	100506658
Components:	KN206468G1 , Occludin gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN206468G2 , Occludin gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN206468BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_001205254](#), [NM_001205255](#), [NM_002538](#)

UniProt ID: [Q16625](#)

Synonyms: BLCPMG; PPP1R115

Summary: This gene encodes an integral membrane protein that is required for cytokine-induced regulation of the tight junction paracellular permeability barrier. Mutations in this gene are thought to be a cause of band-like calcification with simplified gyration and polymicrogyria (BLC-PMG), an autosomal recessive neurologic disorder that is also known as pseudo-TORCH syndrome. Alternative splicing results in multiple transcript variants. A related pseudogene is present 1.5 Mb downstream on the q arm of chromosome 5. [provided by RefSeq, Apr 2011]



[View online »](#)

Product images:

