

Product datasheet for **KN206556LP**

Myocilin (MYOC) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	Myocilin
Locus ID:	4653
Components:	KN206556G1 , Myocilin gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN206556G2 , Myocilin gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN206556LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_000261
UniProt ID:	Q99972
Synonyms:	GLC1A; GPOA; JOAG; JOAG1; myocilin; TIGR
Summary:	MYOC encodes the protein myocilin, which is believed to have a role in cytoskeletal function. MYOC is expressed in many ocular tissues, including the trabecular meshwork, and was revealed to be the trabecular meshwork glucocorticoid-inducible response protein (TIGR). The trabecular meshwork is a specialized eye tissue essential in regulating intraocular pressure, and mutations in MYOC have been identified as the cause of hereditary juvenile-onset open-angle glaucoma. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

