

Product datasheet for **KN206738LP**

KIAA0323 (KHNYN) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	KIAA0323
Locus ID:	23351
Components:	KN206738G1, KIAA0323 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN206738G2, KIAA0323 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN206738LPD, donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003, scramble sequence in pCas-Guide vector
RefSeq:	NM_001290256 , NM_001290257 , NM_015299
UniProt ID:	O15037
Synonyms:	KIAA0323
Summary:	This gene encodes a protein containing a ribonuclease NYN domain. The function of this protein has yet to be determined. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014]



[View online »](#)

Product images:

