

Product datasheet for **KN206965RB**

B3GNT5 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	B3GNT5
Locus ID:	84002
Components:	KN206965G1 , B3GNT5 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN206965G2 , B3GNT5 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN206965RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_032047
UniProt ID:	Q9BYG0
Synonyms:	B3GN-T5; beta3Gn-T5
Summary:	This gene encodes a member of the beta-1,3-N-acetylglucosaminyltransferase family. This enzyme is a type II membrane protein. It exhibits strong activity to transfer GlcNAc to glycolipid substrates and is identified as the most likely candidate for lactotriaosylceramide synthase. This enzyme is essential for the expression of Lewis X epitopes on glycolipids. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

