

Product datasheet for **KN207345LP**

P2RY5 (LPAR6) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	P2RY5
Locus ID:	10161
Components:	KN207345G1 , P2RY5 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN207345G2 , P2RY5 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN207345LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001162497 , NM_001162498 , NM_005767
UniProt ID:	P43657
Synonyms:	ARWH1; HYPT8; LAH3; P2RY5; P2Y5
Summary:	The protein encoded by this gene belongs to the family of G-protein coupled receptors, that are preferentially activated by adenosine and uridine nucleotides. This gene aligns with an internal intron of the retinoblastoma susceptibility gene in the reverse orientation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2009]



[View online »](#)

Product images:

