

Product datasheet for **KN207406RB**

EXTL2 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	EXTL2
Locus ID:	2135
Components:	KN207406G1 , EXTL2 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN207406G2 , EXTL2 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN207406RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001033025 , NM_001261440 , NM_001261441 , NM_001261442 , NM_001439 , NR_048570
UniProt ID:	Q9UBQ6
Synonyms:	EXTR2
Summary:	Glycosyltransferase required for the biosynthesis of heparan-sulfate and responsible for the alternating addition of beta-1-4-linked glucuronic acid (GlcA) and alpha-1-4-linked N-acetylglucosamine (GlcNAc) units to nascent heparan sulfate chains.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

