

## Product datasheet for **KN207967LP**

### **N WASP (WASL) Human Gene Knockout Kit (CRISPR)**

#### **Product data:**

<b>Product Type:</b>	Knockout Kits (CRISPR)
<b>Format:</b>	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
<b>Donor DNA:</b>	Luciferase-Puro
<b>Symbol:</b>	N WASP
<b>Locus ID:</b>	8976
<b>Components:</b>	<b>KN207967G1</b> , N WASP gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN207967G2</b> , N WASP gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN207967LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
<b>Disclaimer:</b>	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
<b>RefSeq:</b>	<a href="#">NM_003941</a>
<b>UniProt ID:</b>	<a href="#">O00401</a>
<b>Synonyms:</b>	N-WASP; NWASP; WASPB
<b>Summary:</b>	This gene encodes a member of the Wiskott-Aldrich syndrome (WAS) protein family. Wiskott-Aldrich syndrome proteins share similar domain structure, and associate with a variety of signaling molecules to alter the actin cytoskeleton. The encoded protein is highly expressed in neural tissues, and interacts with several proteins involved in cytoskeletal organization, including cell division control protein 42 (CDC42) and the actin-related protein-2/3 (ARP2/3) complex. The encoded protein may be involved in the formation of long actin microspikes, and in neurite extension. [provided by RefSeq, Jul 2013]



[View online »](#)

Product images:

