

Product datasheet for **KN208007LP**

NADPH oxidase 4 (NOX4) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	NADPH oxidase 4
Locus ID:	50507
Components:	KN208007G1 , NADPH oxidase 4 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN208007G2 , NADPH oxidase 4 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN208007LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001143836 , NM_001143837 , NM_001291926 , NM_001291927 , NM_001291929 , NM_001300995 , NM_016931 , NR_026571 , NR_120406
UniProt ID:	Q9NPH5
Synonyms:	KOX; KOX-1; RENOX
Summary:	This gene encodes a member of the NOX-family of enzymes that functions as the catalytic subunit the NADPH oxidase complex. The encoded protein is localized to non-phagocytic cells where it acts as an oxygen sensor and catalyzes the reduction of molecular oxygen to various reactive oxygen species (ROS). The ROS generated by this protein have been implicated in numerous biological functions including signal transduction, cell differentiation and tumor cell growth. A pseudogene has been identified on the other arm of chromosome 11. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jan 2009]



[View online »](#)

Product images:

