

Product datasheet for **KN208030RB**

BLAP75 (RMI1) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	BLAP75
Locus ID:	80010
Components:	KN208030G1 , BLAP75 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN208030G2 , BLAP75 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN208030RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_024945 , NM_001358291 , NM_001358292 , NM_001358293 , NM_001358294
UniProt ID:	Q9H9A7
Synonyms:	BLAP75; C9orf76; FAAP75
Summary:	RMI1 is a component of protein complexes that limit DNA crossover formation via the dissolution of double Holliday junctions (Raynard et al., 2006 [PubMed 16595695]).[supplied by OMIM, Mar 2008]



[View online »](#)

Product images:

