

Product datasheet for **KN208556BN**

Heparanase 1 (HPSE) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Heparanase 1
Locus ID:	10855
Components:	KN208556G1 , Heparanase 1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN208556G2 , Heparanase 1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN208556BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001098540 , NM_001166498 , NM_001199830 , NM_006665
UniProt ID:	Q9Y251
Synonyms:	HPA; HPA1; HPR1; HPSE1; HSE1
Summary:	Heparan sulfate proteoglycans are major components of the basement membrane and extracellular matrix. The protein encoded by this gene is an enzyme that cleaves heparan sulfate proteoglycans to permit cell movement through remodeling of the extracellular matrix. In addition, this cleavage can release bioactive molecules from the extracellular matrix. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]



[View online »](#)

Product images:

