

Product datasheet for **KN208818BN**

BCDIN3D Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	BCDIN3D
Locus ID:	144233
Components:	KN208818G1 , BCDIN3D gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN208818G2 , BCDIN3D gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN208818BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_181708](#)

UniProt ID: [Q7Z5W3](#)

Synonyms: BCDIN3 domain containing

Summary: This gene encodes an RNA methyltransferase which belongs to the rossmann fold methyltransferase family, and serves as a 5'-methylphosphate capping enzyme that is specific for cytoplasmic histidyl tRNA. The encoded protein contains an S-adenosylmethionine binding domain and uses the methyl group donor, S-adenosylmethionine. This gene is overexpressed in breast cancer cells, and is related to the tumorigenic phenotype and poor prognosis of breast cancer. The encoded protein is thought to promote the cellular invasion of breast cancer cells, by downregulating the expression of tumor suppressor miRNAs through the dimethylation of the 5-monophosphate of the corresponding precursor miRNAs. [provided by RefSeq, Apr 2017]



[View online »](#)

Product images:

